

Labeling and Hybridization of genomic DNA to NKI cDNA arrays

Maartje Vogel version 26/05/2004

Please check <http://www.nki.nl/nkidep/vansteensel> for updated versions of this protocol.

Measure concentration

Measure the concentration of the cleaned DamID-PCR products on the Fluorimeter using herring sperm as a standard.

DNA sample	Final DNA concentration (ng/ul)

Labeling for microarray hybridization

Genomic DNA can be labeled with a simple random-priming protocol based on Gibco/BRL's Bioprime DNA Labeling kit. This kit is a convenient and inexpensive source of random octamers, reaction buffer, and high concentration Klenow. Do not use the dNTP mix provided in the kit; prepare yourself using PCR grade dNTPs.

10X dNTP mix (store at -20°C):

1.2 mM each dATP, dGTP, and dTTP

0.6 mM dCTP

T₁₀E₁ (10 mM Tris pH 7.5, 1mM EDTA, filtered)

eppie	DNA name	... µl (=...µg)	+ ...µl ddH ₂ O	...µl Cy5/Cy3
1				
2				
3				
4				
5				
6				
7				
8				

- Add 2 µg DNA of the sample to be labeled to a PCR tube.
(Note: for For DamID with lentiviruses add 1 ug.)
- Add ddH₂O to bring the total volume to 21 µl.

- Then add 20 µl of 2.5X Random Primers Solution (from the Bioprime kit Lot no:.....)
- (OR: make a mastermix of DNA, ddH₂O and primer solution)
- Boil 5 min (in PCR machine), and place on ice.
- On ice, add 5 µl 10X dNTP mix.
- Add 3 µl Cy5-dCTP or Cy3-dCTP (Amersham, 1 mM stocks).

Cy5 dCTP Lot no	Cy3 dCTP Lot no
------------------------	------------------------

- Add 1 µl Klenow Fragment (40 U/µl from the Bioprime kit, Lot no:....., Exp. Date:.....).
- Incubate 1.5 hours at 37°C in PCR machine.
- Stop the reaction by adding 5 µl 0.5 M EDTA pH 8.0 (stopbuffer from the Bioprime kit).

Pre-hybridisation (part 1)

Pre-hybridisation buffer (5x SSC, 0.1x SDS and 1% BSA)

(Make a 10% BSA stock, filter and store at –20°C, thaw only once)

For 30 ml:	7.5 ml	20x SSC (filtered)	(60 ml:	15 ml
	150 µl	20% SDS (filtered)		300 µl
	3 ml	10% BSA		6 ml
	19.35 ml	ddH ₂ O		38.7 ml)

- Warm pre-hyb buffer to 42°C prior to use.
- Add the buffer to the 'pre-hyb chamber' (plastic box) and place the slide in the chamber, making sure the array is submerged.
- Pre-hybridize the microarray at 42°C for at least 1 hour (up to 3 hours).

Probe cleanup

Make sure you have:

- TE = T₁₀E₁ (10mM TrisCl pH7.5, 1mM EDTA) filtered.
- Heat block 100 °C (preferably also a 42°C heat block).
- Water bath at 42°C.

- Pre-spin chromaspin column for 3 minutes at 700g (in eppendorf 5417) (see manufactures protocol for instructions). Immediately continue with the next step do not let the column dry out!
- Place a fresh tube below the chromaspin column. Pool the Cy3 and Cy5 reactions and load the sample to the top of the column, make sure that you pipette the probe slowly in the center of the column material, which may come loose from the inside of the column because of spinning.
- Spin column for 3 min at 700 g.

- Collect tube and mark it properly, keep it as dark as possible.
- Mark a microcon column and add 450 µl TE to the column. Pre-spin 5 min at 13.000 g (11.063 rpm) to wash column, discard flow trough. Check if column surface is still wet (if not use new column).
- Apply 200 µl TE to Microcon column.
- Add 64 µl of blocking solution to the Microcon column
- Optional: Add 25 µg pre-digested competitor plasmid to the Microcon column. Not necessary for DamID with lentivirus

Comp. plasmid (..... µg/µl) add µl per hybridization
--

- Load the eluate from the chromaspin column to the Microcon column.

Blocking solution (store at -20°C):

4 µl	20 µg poly d(A) (5 mg/ml Roche)
20 µl	100 µg yeast t-RNA (5 mg/ml stock Life)
40 µl	40 µg human Cot-1 DNA (1 mg/ml Invitrogen)

- Spin Microcon column for 5 min at 13.000g (11.063 rpm) (discard flow trough).
- Add 300 µl TE to top of column and spin 4 min at 13.000 g (discard flow trough).
- Add another 300 µl TE to top of column and spin 5 min at 13.000 g.
- Spin to reduce the volume to 20 µl (judge by eye) check every minute. (When half of the column surface is wet, the volume is <10 µl. (Inspect the color of the labeled DNA by eye).
- Collect probe by inverting the column in a fresh tube (cut of the lid and keep it to seal afterwards) and spin 2 min at 13.000 g.
- Measure the amount of probe and adjust to 20 µl with TE.

Eppie (array number)	
µl probe +µl TE = 20µl
µl probe +µl TE = 20µl
µl probe +µl TE = 20µl
µl probe +µl TE = 20µl
µl probe +µl TE = 20µl
µl probe +µl TE = 20µl

- Heat for 2 minute at 100 °C (dry heat block) and snap cool on ice for 2 minutes. Store the probe for at least 5 minutes at 42°C.
- Make fresh 2x F-hybridization buffer, heat the buffer to 42°C for at least 5 minutes and add 20 µl to the probe, mix the 2x F-buffer well before adding. Keep at 42°C.

2x F-hybridization buffer (pipet in the following order):

500 µl 100% formamide
500 µl 20x SSC
10 µl 20% SDS

1x F-hyb buffer:

dilute 2x F-buffer with ddH₂O

- Place a clean lifter slip on an untreated object glass slide. Determine the volume under the lifter slip with ddH₂O. Adjust the probe volume to the volume under the lifter slip with 1xF-buffer.

Pre-hybridisation (part 2)

- Take the arrays from the pre-hyb buffer, dip them in ddH₂O and collect them in a glass color box with ddH₂O.
- Wash of the pre-hyb buffer for 10 mins in the glass color box.
- Transfer the slide to a second distilled water bath for another 5 mins.
- Then transfer the slide to 100% iso-propanol for 5 minutes in a third glass color box.
- (Air dry or) briefly spin dry the arrays for ~2 mins at 600 rpm. Make sure that the time between the drying and the actual hybridization is not longer than 15 minutes.

Hybridization

- Clean the lifter slip with soap, plenty of water, ethanol and acetone. Dirty lifter slips give high background. (When the slips are not used before: start to wash with acetone.)
- Place the array in the chamber with the array site up; make sure that the array does not touch the inner sides of the chamber.
- Put the lifter slip on the array with the with Teflon linings facing down, so the lifter slip is not directly touching the glass.

The next 3 steps should be performed as quickly as possible to prevent evaporation of the probe.

- Spin probe for 1min at max speed to get rid of floating precipitates.
- Apply ~40 ul 1x F-hybridization buffer to the far ends of the chamber.
- Adjust probe to array, the probe will be sucked under the cover slide, due to capillary force. (Inject the probe slowly from one corner of the lifter slip until array is thinly covered – the pipette tip should be as close to the slip as possible without touching it)
- Close the hyb chamber, keep it horizontally, move it only slowly and make sure that the array does not move in the chamber.
- Incubate at 42°C (in a stove) for 4x o/n.

Washes

Prepare the washing solutions, the first three washing solutions are kept at 42°C, solution 4 and 5 are kept at RT.

WASH 1: 5x SSC / 0.1%SDS (149 ml H₂O +50 ml 20xSSC + 1 ml 20% SDS)
WASH 2: 2x SSC / 0.1%SDS (179 ml H₂O +20 ml 20xSSC + 1 ml 20% SDS)
WASH 3: 1x SSC (190 ml H₂O +10 ml 20xSSC)
WASH 4: 0.2x SSC (198 ml H₂O +2 ml 20xSSC)
WASH 5: 0.05x SSC (200 ml H₂O +0.5 ml 20xSSC)

- Take the array from hyb chamber and submerge it as quickly as possible in the first wash solution. (Process – put washing buffer 1 in a 50 ml falcon tube, submerge array standing upright, lifter slip will slide off after a while without scratching the array.)
- Wash the array in solution 2 for about half a minute, plunging gently.
- Wash the array in solution 3 for about 5 mins, also plunging gently. (Arrays can be collected in a rack in solution 3).
- Wash the array(s) in solution 4 for about 2 mins, again plunging gently.
- Wash array(s) in solution 5 for 20 seconds without plunging and spin dry quickly (5 minutes 600 rpm). Be careful to use a good counter weight. DO NOT LET SLIDES AIR DRY. Keep the dry arrays dark.
- Scan the array (arrays can be stored for three to four days in the dark before scanning).

General notes

Concentration measurements are preferably done with a Fluorimeter (DyNA 200 from Hoefer) with Hoechst using herring sperm DNA for standard calibration. The DNA amounts used in this protocol are based on Hoechst measurements using herring sperm standards.

Concentration measurements can be done on the nanodrop spectrophotometer. Use half of the amount of DNA indicated in the protocol when the nanodrop is used.
For example: Dpn1 digest 2.5 ug gDNA (nanodrop) instead of 5 ug (Hoechst).

Ordering information

Labeling and microarray hybridization:

Invitrogen	BioPrime DNA labeling System	30 rxns	18094-011
Invitrogen	yeast tRNA	25 mg	15401-011
Amersham	Cy3-dCTP	25 nmol	PA53021
Amersham	Cy5-dCTP	25 nmol	PA55021
Roche	Poly [A]	100mg	108626
Millipore	Microcon YM-30 filters	100	42410
Invitrogen	human Cot-1 DNA	500ug	15279-011
Clontech	Cromaspin +TE-30 columns	50	636069