




www.roche-applied-science.com



Genome Sequencer 20 System

Accelerate Amplicon Analysis



Genome Sequencer 20 Amplicon Software screen.
The histogram shows the deviations between read consensus sequences and reference sequences. In addition, the software displays the corresponding alignment of reads and reference sequences for comparison purposes.

Visit www.genome-sequencing.com for more information and published references, call **800 262 4911** (U.S.), or contact your local sales representative.

Clonal sequencing of amplicons enables ultra-deep genetic analysis at levels of sensitivity not previously possible.

Cancer Research

- Discover somatic mutations with unprecedented speed, accuracy, and sensitivity.*
- Study DNA-methylation patterns with high precision.

Medical Resequencing

- Uncover and screen SNPs at the population and individual level.
- Deconvolute complex populations of sequence data.
- Eliminate tedious and expensive manual analysis of mixed chromatograms.

First to the finish with ultra-deep analysis of amplicons – new for the Genome Sequencer 20 System.

*Thomas, R. *et al.*, Nature Medicine advance online publication, 25 June 2006 (doi:10.1038/nm1437)

454 LIFE SCIENCES



Diagnostics

THE CONDENSED PROTOCOLS

FROM MOLECULAR CLONING: A LABORATORY MANUAL

By Joseph Sambrook, *Peter MacCallum Cancer Institute, Melbourne, Australia*, and David Russell, *University of Texas Southwestern Medical Center, Dallas*

The Condensed Protocols From Molecular Cloning: A Laboratory Manual is a single-volume adaptation of the three-volume third edition of *Molecular Cloning: A Laboratory Manual*. This condensed book contains only the step-by-step portions of the protocols, accompanied by selected appendices from the world's best-selling manual of molecular biology techniques. Each protocol is cross-referenced to the appropriate pages in the original manual. This affordable companion volume, designed for bench use, offers individual investigators the opportunity to have their own personal collection of short protocols from the essential *Molecular Cloning*.


2006, 800 pp., illus., appendices, index

Hardcover \$150

Paperback \$95

ISBN 0-87969-772-5

ISBN 0-87969-771-7



REVIEWS OF *MOLECULAR CLONING: A LABORATORY MANUAL*, 3rd ed.

“Any basic research laboratory using molecular biology techniques will benefit from having a copy [of *Molecular Cloning*] on hand. . . Originally published almost 20 years ago, the first two editions of this book have been staples of molecular biology with a proven reputation for accuracy and thoroughness.”
—*The Scientist*

“The next generation of *Molecular Cloning* not only carries on the proud heritage of the first two editions but also admirably expands on that tradition to provide a truly essential laboratory manual.”
—*Trends in Microbiology*

“This book should be the first port of call when establishing a molecular technique new to a laboratory.”
—*Trends in Cell Biology*

“. . . with this new edition and its companion website, *Molecular Cloning* has once again established its primacy as the molecular laboratory manual and is likely to be found on lab benches, and book-marked on web browsers, around the world.”
—*Trends in Neurosciences*

Contents

- | | | | |
|--|---|---|--|
| 1. Using Plasmid Vectors in Molecular Cloning | 6. Preparation and Analysis of Eukaryotic Genomic DNA | 10. Working with Synthetic Oligonucleotide Probes | 16. Introducing Cloned Genes into Cultured Mammalian Cells |
| 2. Bacteriophage λ and Its Vectors | 7. Extraction, Purification, and Analysis of mRNA from Eukaryotic Cells | 11. Preparation of cDNA Libraries and Gene Identification | 17. Analysis of Gene Expression in Mammalian Cells |
| 3. Working with Bacteriophage M13 Vectors | 8. In Vitro Amplification of DNA by the Polymerase Chain Reaction | 12. DNA Sequencing | 18. Protein Interaction Technologies |
| 4. Working with High-capacity Vectors | 9. Preparation of Radiolabeled DNA and RNA Probes | 13. Mutagenesis | Appendices |
| 5. Gel Electrophoresis of DNA and Pulsed-field Agarose Gel Electrophoresis | | 14. Screening Expression Libraries | Index |
| | | 15. Expression of Cloned Genes in <i>Escherichia coli</i> | |

To order or request additional information:

Call: 1-800-843-4388 (Continental US and Canada) 516-422-4100 (All other locations)

FAX: 516-422-4097

E-mail: cshpress@cshl.edu or WWW Site <http://www.cshlpress.com/>

Write: Cold Spring Harbor Laboratory Press, 500 Sunnyside Blvd, Woodbury, NY 11797-2924





focus on discovery.

Discover More with the World's Finest Full-Length Clone Collection

The Mammalian Gene Collection (MGC) is quite simply the best clone collection in the world, brought to you by the greatest minds in genomic research. Generated and sequenced by a consortium of leading NIH funded laboratories, our cDNA resource contains a complete open reading frame (ORF) for every human, mouse, and rat gene represented. The power of discovery is at your fingertips.

- Every MGC clone is fully sequenced with complete ORF
- Assay-ready DNA in microplate format upon request
- Online clone query quickly identifies the clone you need

"Open Biosystems' MGC collection is an invaluable resource for CCSB, the Dana-Farber Center for Cancer Systems Biology. It has cut years off the development of a human protein-protein interaction map."

Marc Vidal
Dana-Farber Cancer Institute

Why experiment when you can discover?

Your expert partner for
Genomic, RNAi and Proteomic resources.



Toll Free 888.412.2225
www.openbiosystems.com

GENOMICS • RNAi • PROTEOMICS
express repress detect



yield of dreams.


Taq DNA Polymerase from New England Biolabs

HIGH YIELD, ROBUST AND RELIABLE PCR REACTIONS IN CONVENIENT FORMATS


Looking for the right solution for your high yield PCR? Choose recombinant *Taq* DNA Polymerase from New England Biolabs. As the leader in enzyme technology, New England Biolabs provides the highest quality recombinant *Taq* at exceptional value. Our expanded selection of *Taq* based products includes kits, master mixes, and a choice of reaction buffers. Choose *Taq* DNA polymerase from NEB for guaranteed Performance, Convenience and Results.

- **Taq with Standard Buffer**  **M0273S/L**
Compatible with existing assay systems and high throughput applications, is always detergent free and available with or without $MgCl_2$
- **Taq with ThermoPol Buffer**  **M0267S/L**
Promotes high yields under demanding conditions, available with or without $MgSO_4$ and detergent
- **Taq 2X Master Mix**  **M0270S/L**
Just add template and primers
- **NEW Quick-Load™**
Taq 2X Master Mix  **M0271S/L**
Load PCR products directly onto agarose gels
- **Taq PCR Kits**  **E5000S/E5100S**
Contains reagents for 200 PCR reactions, available with or without controls

 = Recombinant



Selective and specific PCR amplification from Human Genomic DNA:
Specific amplicons of 0.5, 1, 2, and 4 kb from human genomic DNA were amplified by *Taq* DNA Polymerase with Standard Buffer for 30 cycles. Marker (M) shown is 2-Log DNA Ladder (NEB #N3200).



Versatility of the Taq 2X Master Mix:
30 ng human genomic DNA (hDNA) or 0.1 ng lambda DNA (λ DNA) was amplified in the presence of 200 nM primers in a 25 μ l volume. Marker (M) shown is 2-Log DNA Ladder (NEB #N3200).

For more information and international distribution network, please visit www.neb.com

- **New England Biolabs Inc.** 240 County Road, Ipswich, MA 01938 USA 1-800-NEB-LABS Tel. (978) 927-5054 Fax (978) 921-1350 info@neb.com
- **Canada** Tel. (800) 387-1095 info@ca.neb.com
- **Germany** Tel. 0800/246 5227 info@de.neb.com
- **UK** Tel. (0800) 318486 info@uk.neb.com
- **China** Tel. 010-82378266 beijing@neb-china.com

 NEW ENGLAND
BioLabs® Inc.
the leader in enzyme technology

LESS
WASTE.

MORE
BAND.



PCR so reliable, we guarantee your results. Get successful PCR the first time with AccuPrime™ enzymes and primers from Invitrogen. Put an end to waste: reduce rework, optimization, and repeat reactions. Start getting used to PCR you can count on. And if something goes wrong, we'll troubleshoot it and make it right—even redesign and synthesize primers at no charge. See how dependable PCR can be at www.invitrogen.com/accuprime.



invitrogen™

BIOSUPPLYNET

A Current, Integrated Information Source
for Life Science Laboratory Supplies



BIOSUPPLYNET

Find A Kit

Many suppliers offer time and effort saving kits to perform research techniques. Search our exclusive database of available research kits.

BioToolKit

An annotated directory of over 1000 links to online molecular biology resources, including basic research tools and advanced applications for genome, transcriptome, and proteome data retrieval, analysis, and visualization.

Featured Subject Areas

Download free protocols, background information and find suppliers for reagents and equipment for the hottest fast-growing research methods and areas of interest.

Featured Protocols

Download free protocols from Cold Spring Harbor Laboratory Press manuals, and find sources for the reagents and equipment needed.

Career Center

Find or place listings for available positions in academic research, pharmaceutical companies, clinical research facilities, and the biotechnology industry.

Searching

Searching is very convenient. Info on over 20,000 products from over 6,500 suppliers is right at your finger tips.


CSHL Press Announces the New and Improved BioSupplyNet.com!

- CSHL Press has relaunched BioSupplyNet, adding new features and improved functionality
- SAVE TIME! No need to search through stacks of catalogs and numerous websites
- Over 20,000 products and supplies listed (indexed by type and/or by company)
- Over 6,500 prominent suppliers participating and the number continues to grow
- Searching is FAST and EASY
- NO REGISTRATION necessary

In addition to searching for products, Scientists can quickly:

- Get immediate inside info about new products and special deals
- Find kits to perform research techniques
- Visit Featured Subject Areas containing background info, free protocols, hottest fast-growing research methods and more! Now featuring PCR and RNAi
- Download free protocols from Cold Spring Harbor Laboratory Press manuals
- Visit our Career Center to search and post job listings
- Get up to the minute Science news
- Order Catalogs
- Sign up for Newsletters

Visit BioSupplyNet Today!



www.biosupplynet.com

CSH Protocols Is Now Online!



- A new, online source of trusted techniques in molecular and cellular biology
- Contains new and classic protocols presented step-by-step with cautions and troubleshooting
- Frequently updated and annotated
- Interactive, customizable, and fully searchable

Cold Spring Harbor Laboratory is renowned for its teaching of biomedical research techniques. For decades, participants in its celebrated, hands-on courses and users of its laboratory manuals have gained access to the most authoritative and reliable methods in molecular and cellular biology. Now that access has moved online.

Visit *Cold Spring Harbor Protocols* today and discover a rich, interactive source of new and classic research techniques. The site is fully searchable, with many tools that can be customized by users including topic-based alerting and personal folders. Through a web-based editorial process, users also have the opportunity to add refereed comments to each protocol. Links in the online protocols offer additional resources, and step-by-step instructions print out in a convenient form, complete with materials, cautions, and troubleshooting advice. Each protocol is citable and presented and edited in the style that has made *Molecular Cloning*, *Antibodies*, *Cells*, and many other CSH manuals essential to the work of scientists worldwide. The collection of more than 900 protocols released in 2006 will be continuously expanded, updated, and annotated by the originators and users of the techniques.

CSH Protocols is created by Cold Spring Harbor Laboratory Press in association with HighWire Press of Stanford University.

Cold Spring Harbor Protocols • Executive Editor: Dr. Michael Ronemus
ISSN 1559-6095 / online, monthly • Available exclusively via institutional site license

Subject Coverage

Antibodies	Genetics	Newly Added Protocols
Bioinformatics/Genomics	High Throughput Analysis	Plant Biology
Cell Biology	Imaging/Microscopy	Polymerase Chain Reaction (PCR)
Chromatography	Immunology	Proteins and Proteomics
Computational Biology	Laboratory Organisms	RNA Interference (RNAi)/siRNA
DNA Delivery/Gene Transfer	Molecular Biology	Stem Cells
Electrophoresis	Neuroscience	Transgenic Technology

*CSH Protocols—continuing Cold Spring Harbor Laboratory's
60-year tradition as a source of trusted techniques.*

Visit Today!
www.cshprotocols.org

For pricing information or to request a free trial:

Call: 1-800-843-4388 (Continental US and Canada) 516-422-4100 (All other locations)

FAX: 516-422-4097

E-mail: cshpress@cshl.edu

Write: Cold Spring Harbor Laboratory Press, 500 Sunnyside Blvd, Woodbury, NY 11797-2924






Get *Dependable Delivery* of
 High Sensitivity, High Yield, High Fidelity & Reproducible
RNA Amplification

Whether you use GeneChip® or other arrays, EPICENTRE Biotechnologies' **TargetAmp™ RNA Amplification Kits** deliver dependable, high quality microarray results.


High Yield

Input Total RNA	aRNA Yield	Fold-amplification
TargetAmp™ 1-Round aRNA Amplification Kit 103		
25 ng	2.9 µg	5800
100 ng	11.4 µg	5700
500 ng	69.9 µg	6990
TargetAmp™ 2-Round aRNA Amplification Kit 2.0		
1 cell (10 pg)	3 µg	6 x 10 ⁶
10 cells (100 pg)	53 µg	10.6 x 10 ⁶
50 cells (500 pg)	161 µg	6.4 x 10 ⁶

High Fidelity




Reproducible



High Quality Microarray Results

	TargetAmp™ 1-Round Kits	TargetAmp™ 2-Round Kits
Fold-Amplification	>5,000	>5,000,000
3'/5' ratio GAPDH	1.05 ±0.05	5.6 ±1.5
3'/5' ratio β-actin	2.8 ±0.01	10.7 ±0.07
Present Calls	56.3%	50.7%



- TargetAmp™ 1-Round Aminoallyl-aRNA Amplification Kit 101
- TargetAmp™ 1-Round aRNA Amplification Kit 103
- TargetAmp™ 2-Round Aminoallyl-aRNA Amplification Kit 1.0
- TargetAmp™ 2-Round aRNA Amplification Kit 2.0

For more information, go to www.EpiBio.com and enter this QuickInfo code: **TZAX4**



This product is covered by intellectual property licensed to EPICENTRE Technologies Corporation from Johnson & Johnson Pharmaceutical Research & Development, L.L.C. This product is covered by U.S. Patents licensed exclusively to Incyte Corporation and sublicensed to EPICENTRE Technologies Corporation. See www.EpiBio.com website for complete license statements. GeneChip is a registered trademark of Affymetrix. TargetAmp is a trademark of EPICENTRE.

CAREER TRACKS

Dedicated entirely to Employment, Conferences, Meetings, Fellowships, and Grants



Register by September 15th and Save up to \$200!

 Cambridge Healthtech Institute's *Fourth Annual*

Discovery on Target 2006

Target Discovery & Target Validation Strategies

October 23-26, 2006 • World Trade Center • Boston, MA



"Highly focused, very insightful, interactive and thought provoking meeting."
Venata Sabbiseth, Brigham and Women's Hospital

CONCURRENT CONFERENCES

PRECLINICAL DISEASE MODELS

To Identify and Validate Novel Targets

RNAi

From Target Discovery & Validation to Therapeutic Development

EXECUTIVES ON TARGET

Maximizing the Impact of Drug Discovery

CHEMICAL GENOMICS

Small Molecules to Discover Targets

ION CHANNELS

New

An Emerging Target

PLUS! PRE-CONFERENCE EVENTS

Signature Sponsor:  BANNER & WITCOFF, LTD.
INTELLECTUAL PROPERTY LAW

Corporate Sponsors:



VASTOX

Reference Keycode S54AD when registering online

www.DiscoveryOnTarget.com



Cambridge Healthtech Institute 250 First Avenue, Suite 300, Needham, Massachusetts 02494 • T: 617-630-1300 or toll-free in the U.S. 888-999-6288 • F: 617-630-1325 • www.healthtech.com



Plant & Animal Genome XV
The International Conference on the
Status of
Plant & Animal Genome Research
January 13-17, 2007
Town & Country Hotel
San Diego, California

Organizing Committee

Chairman:

Stephen R. Heller, NIST, USA

Plants:

Michael Gale, John Innes Center, UK
Ed Kaleikau, USDA/CSREES, USA
Dave Matthews, USDA, ARS Cornell University, USA
Graham Moore, John Innes Center, UK
Jerome P. Miksche, Emeritus Director, USDA Plant
Genome Program, USA
Rod Wing, University of Arizona, USA

Animals:

Mary Delany, USA
Michel Georges, University of Liege, Belgium
Clare A. Gill, Texas A&M University, USA
Jim Reecy, USA

Abstract Coordinators:

Victoria Carollo, USDA, ARS, WRRRC, USA
Gerard Lazo, USDA/ARS/WRRRC, Albany, CA, USA
David Grant, USDA/ARS & Iowa State University, USA

Sponsors

USDA, Agricultural Research Service
USDA, National Agricultural Library
USDA, NRI Competitive Grants Office
USDA, Cooperative State Research, Education, and
Extension Service (CSREES)
Rockefeller Foundation
John Innes Centre
NCGR, National Center for Genome Resources

Speakers

Martina Newell-McGloughlin,
University of California - Davis
Ian T. Baldwin,
Max Planck Institute for Chemical Ecology - Jena
Peter J. Hunter, University of Auckland
David Bartel, MIT
Susan McCouch, Cornell University
Detlef Weigel, Max Plank Institute - Tübingen
Kirsten Lindblad-Toh, Harvard/MIT Broad Institute

Workshops

Abiotic Stress / Allele Mining / Apomixis Aquaculture / Banana
(Musa) Genomics / Barley / Bioinformatics / Brassicas /
Brachypodium Distachyon / Cattle/Sheep / Challenge Program:
Unlocking Crop Genetic Diversity for the Poor / Citrus
Compositae / Computer Demonstrations / Cotton / Compositae /
Connectrons Cool Season Legumes / Curcubit / Equine / Forage &
Turf Plants / Forest Trees Fruit and Nut Crops / Functional Genom-
ics / Host Pathogen Interactions / Insect Genetics / ICSB / ICGI /
Int'l Grape Genome Project / Int'l Lolium Genome
Initiative / IGGI / ITMI / Large-Insert DNA Libraries and Their
Applications / Legumes/ Maize / Microarray Analysis /
Molecular Markers for Plant Breeders/ Mutation Screening /
NRSP-8 / NC 1010 Development and Implementation of
Ontologies in the Database / Organellar Genetics / Plant
Cytogenetics / Plant Development and Signal Networks /
Interagency Working Group on Plant Genomics / Plant
Interactions with Pests and Pathogens / Plant Transgene
Genetics / Plant Reproductive Genomics / Polyploidy / Poultry /
Proteomics / QTL Cloning / Reduced-representation / Sequencing
Methods and Applications / Rice / Rice Blast /
Root Genomics / Solanaceae / Sorghum and Millets / Soybean Ge-
nomics / Statistical Genomics / Sugar Beet / Swine /
Swine Genome Sequencing / TAIR /
Weedy and Invasive Plant Genomics

Organizer

Scherago International
525 Washington Blvd., Ste. 3310
Jersey City, NJ 07310
201-653-4777 x20
201-653-5705
E-mail: pag@schicago.com

For complete details, including on-line registration, visit our website at
www.intl-pag.org