

JUST RIGHT FOR GENOMICS

ROOTED IN GENOMIC RESEARCH, DEVELOPED FOR THE FUTURE.

By simultaneously quantifying and qualifying nucleic acid samples in parallel, the Fragment Analyzer™ is transforming sample prep analysis for the world's leading genomic research institutions. Automate genomic QC for an array of applications, including total and degraded RNA isolations, genomic DNA extractions and NGS library preparations—giving you better results in less time, with less effort.



More at AATI-US.COM

Expand Your Epigenetics Discovery by Sequencing More of the Methylome

Accel-NGS® Methyl-Seq Library Kit

- Simple 2-hour protocol
- DNA inputs from 100 pg to 100 ng
- Low bias library preparation
- Compatible with WGBS, RRBS, and other enrichment methods



© 2016, Swift Biosciences, Inc. The Swift logo is a trademark and Accel-NGS is a registered trademark of Swift Biosciences. 16-0963, 07/16

The S is for Simplicity

The new Ion S5™ System.
Targeted sequencing has
never been simpler.

Simple library prep tools, cartridge-based reagents and automated data analysis have reduced DNA-to-data hands-on time to less than 45 minutes. So you'll spend less time doing routine molecular biology, and more time informing time-sensitive decisions.



Ion AmpliSeq™ technology

As little as 1 ng low-quality DNA
sample input for library prep



Cartridge-based reagents

Less than 15 minutes of
sequencing setup time




2.5 to 4 hours of run time

Fastest run time of any
benchtop sequencer

Watch the Ion S5 System in action at
thermofisher.com/ionS5

ThermoFisher
SCIENTIFIC



Be part of the CRISPR revolution!

Talk to us about your ideas at:
info@silence-therapeutics.com


DESIGN - ENGINEER - INNOVATE

CRISPR

Your complete line of CRISPR genome editing tools

- The most comprehensive commercial screening portfolio from MilliporeSigma, the leader in lentiviral manufacturing
 - In partnership with Sanger Wellcome Trust, MilliporeSigma offers the first and only whole genome arrayed CRISPR library available for both human and mouse. Streamline your screening workflow with one optimized clone per well
 - GeCKO v2 whole genome CRISPR pools. Maximize statistical confidence with full genome coverage
 - Custom library panels and pools. Create an exclusive Sanger clone array or GeCKO pool tailored to your research
- SygRNA™, synthetic crRNA and tracrRNA, in conjunction with Cas9 protein offer fast turnaround time to accelerate your experimental results
- Advanced technical support to enhance your gene editing efforts. Contact CRISPR@sial.com with your genome editing questions

For more information, visit
sigma-aldrich.com/CRISPR



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

Copyright © 2016 EMD Millipore Corporation. All Rights Reserved. MilliporeSigma is a trademark of Merck KGaA, Darmstadt, Germany and its affiliates.




Next Generation Sequencing (NGS) Services

SBS Genetech provides a wide variety of next generation sequencing suitable for most research needs. We offer comprehensive services including library preparation, sample validation, sequencing, and bioinformatics.

●Comprehensive NGS Services

Genomic DNA Sequencing
Strand Specific RNA Sequencing
miRNA Sequencing
PCR Amplicon Sequencing
16S Metagenomic Sequencing



●Sequencing Only Services

HiSeq Sequencing
MiSeq Sequencing



●Specialized Sequencing Services

Library Generation
Project Planning and Consultation

Interested? Please send your inquiries to jacklan@sbsbio.com





CRISPR

Dharmacon reagents help you move forward. Faster.

Simplify your gene editing workflow with the Dharmacon™ Edit-R™ CRISPR-Cas9 platform. Custom or predesigned ready-to-use DNA and synthetic RNA leverage a validated algorithm to select highly functional, specific targets for gene knockout – enabling you to quickly assess multiple target sites per gene across one or hundreds of genes.

CRISPR Guide RNAs | CRISPR Screening Libraries | Cas9 Nucleases

crRNAs that have high functionality scores show high editing efficiency



HEK293T-CAG-Cas9 cells were transfected with either high-scoring or low-scoring crRNAs (50 nM crRNA:tracrRNA) using DharmaFECT 1 transfection reagent (0.25 μ L/well) in 96-well format. Gene editing efficiencies were determined using next-generation sequencing. 93% of the top 10 high-scoring crRNAs targeting ten different genes have > 40% indel formation and only 33% of the 10 lowest scoring designs have > 40% indel formation.



Works in difficult to reach places.



In 2015, Ebola ravaged western Africa. Using Recombinase Polymerase Amplification (RPA) from TwistDx, a group of ground-breaking researchers created a solar-powered mobile lab that could begin to **detect Ebola virus in 3 minutes** without a thermocycler.

Read how RPA was used in the Ebola crisis
www.goo.gl/vwEfH9

RPA. It really works.
twistdx.co.uk | +44 (0)1223 496700

Scan the QR code for more information on our products and services



Empowering Research Through Smart Solutions

The right tools for your bioinformatics analysis needs

novoAlign

- The market's leading aligner for accurate and speedy multi-threaded sequence alignment.
- Powerful tool designed for mapping of short reads onto a reference genome.
- Now available on BaseSpace.

novoSort

- The fastest multi-threaded sort/merge tools for BAM files.
- Sort and mark duplicates in a flash.

novoLR package

- novoLRcleaver
- novoLRcorrector
- novoLRpolish
- Genome assembly prep using hybrid technology for mixed short reads and single molecule long reads.
- Long reads hybrid correction using in-house algorithms for read ranking and correction.
- Post-assembly sequence polishing.

novoWorx

A on-site, integrative, customizable workbench that allows users to run an entire pipeline without using command lines interface. The platform mainly utilizes Novocraft's proprietary software; novoAlign for alignment and novoSort for sorting and SAM to BAM conversion. A combination of unique in-house softwares and open source modules to decipher your big data into meaningful results.

novoClinic

A patient-centric NGS targeted sequence analysis platform that provides integrated sample tracking for quality control and compliance. The built-in customizable analysis pipeline and straight-forward data reporting system will ease the burden of data mining and interpretation, allowing clinicians to focus on diagnosis and treatment.



NOVOCRAFT TECHNOLOGIES SDN. BHD.

C-23A-05, 3 Two Square, Section 19, 46300 Petaling Jaya, Selangor, Malaysia

+603 7960 0541

office@novocraft.com

+603 7960 0540

http://www.novocraft.com

Novocraft-Technologies-Sdn-Bhd

TRIZOL[®] IN

Direct-zol[™] RNA Purification


- ✓ No phase separation
- ✓ No precipitation
- ✓ No RNA loss




RNA OUT

7 minutes from TRIZOL[®] to RNA:
Mix sample in TRIZOL[®], TRI Reagent[®], or similar.
Add ethanol & load into the Zymo-Spin[™] column.
Patented technology


High-quality small & large RNA



Non-biased miRNA recovery



Highly cited



Learn more at

www.zymoresearch.com/direct-zol

Samples available at zymoresearch.com/direct-zol or call (949) 679-1190

Product	Cat. No.	Size
Direct-zol [™] RNA MiniPrep Plus	R2070, R2071* R2072, R2073*	50 Preps. (≤100 µg RNA) 200 Preps. (≤100 µg RNA)
Direct-zol [™] RNA MiniPrep	R2050, R2051* R2052, R2053*	50 Preps. (≤50 µg RNA) 200 Preps. (≤50 µg RNA)
Direct-zol [™] RNA MicroPrep	R2060, R2061* R2062, R2063*	50 Preps. (≤10 µg RNA) 200 Preps. (≤10 µg RNA)

96-well spin plate & MagBead formats also available

* = Supplied with TRI Reagent[®]



ZYMO RESEARCH
The Beauty of Science is to Make Things Simple[®]

Tel: (949) 679-1190

www.zymoresearch.com



Synthetic Biology • Gene Expression • qPCR • Sequencing • Single-cell Genomics • Cell-based Assays • HTS and Secondary Screening

Echo® Acoustic LIQUID HANDLING

Reduce Reaction Volume

100-Fold

Empower your research with assay miniaturization

The Echo acoustic droplet ejection (ADE) technology is a tipless, non-contact, low-volume liquid dispensing method that allows researchers to work effortlessly in their labs. This powerful nanoliter-scale technology downscales PCRs and DNA assembly methods from the microliter to the nanoliter level with high efficiency – effectively cutting reagent costs and reducing the volume required 100-fold. With the ability to transfer volumes as low as 2.5 nanoliters, the Echo Acoustic Liquid Handler extends the useful life of a primer library and eliminates the need to dilute high concentration primers, saving storage and primer costs.

- Reduce DNA assembly costs
- Quickly pool DNA oligos or fragments
- Eliminate primer dilutions and errors
- Decreased overall operating costs
- No more pipetting



Echo® 555 LIQUID HANDLER

LABCYTE 
The Future of Science is Sound

© 2016 LABCYTE INC. All rights reserved. Labcyte®, Echo®, and the Labcyte logo are registered trademarks or trademarks of Labcyte Inc., in the U.S. and/or other countries.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

 @LabcyteInc
info-us@labcyte.com

MOVING LIQUIDS *with* SOUND

Extraordinary variation of floral organs in *Phalaenopsis* orchid is due to differential expression of multiple genes in sepal, petal, and labellum.¹

¹ Hsiao, Y. et al. *Gene*. 518, 91-100 (2013).

Differential gene expression creates beauty

Differential gene expression plays a significant role in development of many species including orchids.

Phalaenopsis orchids are important species for development and evolutionary studies. Moreover, they are well recognized for their exceptional beauty. The flowers have several organs (i.e. sepal, petal, and labellum), which are significantly different, however together compose conspicuous and harmonious look (view photo). The unique shape of these flower organs is a result of differential expression of multiple genes involved in their development.

Lexogen is focusing on development of accurate and affordable tools for transcriptome analysis with RNA-Seq. QuantSeq is a dedicated kit for expression profiling. It is an easy protocol for producing highly strand-specific next generation sequencing libraries from the 3' end of polyadenylated RNA. Only one fragment per transcript is produced, directly linking the number of mapping reads to the gene expression values. Restricted length saves sequencing space and allows for high level of multiplexing, enabling cost-efficient and fast RNA-Seq experiment.



Expression Profiling RNA-Seq Library Prep Kit

- Gene expression analysis
- Exact 3' UTR tagging
- From 100 pg total RNA input including low quality RNA and FFPE samples
- Cost-effective sequencing of up to 96 samples / lane
- Ready-to-sequence libraries in 4.5 hours
- Illumina™ and Ion Torrent™ compatible
- Reduced data analysis time
- Custom solutions available

LEXOGEN
Enabling complete transcriptome sequencing

www.lexogen.com




The Transfection Experts

X2 for CRISPR

Optimized DNA, RNA
and RNP Delivery
Using CRISPR/Cas9

Download the white paper:
www.mirusbio.com/CRISPR



T7E1 cleavage efficiency in 293T/17 cells transfected with Cas9 pDNA/gRNA or Cas9 protein/gRNA (RNP) using *TransIT-X2*®.

TransIT-X2® Dynamic Delivery System

Now you can use the same transfection reagent for your knockout, knockdown or expression experiments. The *TransIT-X2*® Dynamic Delivery System takes CRISPR genome editing a step further with an advanced polymeric technology that efficiently delivers plasmid DNA, small RNAs such as siRNA and CRISPR guide RNA or RNP complexes.

The *TransIT-X2*® Dynamic Delivery System delivers CRISPR/Cas components in multiple formats:

- DNA**—deliver plasmid DNA expressing Cas9 or guide RNA
- RNA**—deliver sgRNA or crRNA:tracrRNA
- Protein**—deliver Cas9:gRNA RNP complexes

ADVANCE YOUR TRANSFECTIONS.

Request a FREE SAMPLE of
TransIT-X2® Dynamic Delivery System.

Visit www.mirusbio.com, call 888.530.0801 (U.S. only) or +1.608.441.2852

mirusbio.com

Providing gene delivery expertise since 1995

©2016 All rights reserved Mirus Bio LLC. Mirus Bio and *TransIT-X2* are registered trademarks of Mirus Bio LLC.





BD FACSeq™ Cell Sorter and BD™ Precise Assays

Enhanced transcript analysis of single cells

Seamless workflow from sample to digital data

- Acquire protein and transcript information from every cell analyzed
- BD™ Molecular Index technology removes amplification bias, improving gene expression analysis results
- Easily prepare NGS-ready libraries in one tube, in as fast as five hours
- Quickly select and isolate hundreds to thousands of cells


See more at bdbiosciences.com/go/facsseq

BD Life Sciences, 2350 Qume Drive, San Jose, CA 95131

bdbiosciences.com

Class 1 Laser Product. For Research Use Only. Not for use in diagnostic or therapeutic procedures.
© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.
23-17810-01





2016-2017 SCIENTIFIC CONFERENCES

Presenting the most significant research on cancer etiology, prevention, diagnosis, and treatment

EORTC-NCI-EMA-AACR International Conference on Innovation and Biomarkers in Cancer Drug Development

Conference Co-Chairpersons:
Denis A. Lacombe and John W. Martens
September 8-9, 2016 • Brussels, Belgium

11th Biennial Ovarian Cancer Research Symposium

Conference Co-Chairpersons: *Deborah K. Armstrong, Martin M. Matzuk, Gordon B. Mills, and Saul E. Rivkin*
September 12-13, 2016 • Seattle, WA
Co-presented with the Rivkin Center for Ovarian Cancer

Colorectal Cancer: From Initiation to Outcomes

Conference Co-Chairpersons: *Ernest T. Hawk, Steven H. Itzkowitz, Kenneth W. Kinzler, and Johanna W. Lampe*
September 17-20, 2016 • Tampa, FL

Ninth AACR Conference on The Science of Cancer Health Disparities in Racial/Ethnic Minorities and the Medically Underserved

Conference Co-Chairpersons: *Rick A. Kittles, Folakemi T. Odedina, Jeffrey N. Weitzel, and Jun J. Yang*
September 25-28, 2016 • Fort Lauderdale, FL

CRI-CIMT-EATI-AACR Second International Cancer Immunotherapy Conference: Translating Science into Survival

Conference Co-Chairpersons: *James P. Allison, Philip D. Greenberg, Christoph Huber, and Guido Kroemer*
September 25-28, 2016 • New York, NY

Tumor Immunology and Immunotherapy

Conference Co-Chairpersons: *James P. Allison, Pamela S. Ohashi, Antoni Ribas, and Ton Schumacher*
October 20-23, 2016 • Boston, MA

Translational Control of Cancer: A New Frontier in Cancer Biology and Therapy

Conference Co-Chairpersons: *Jennifer A. Doudna, Frank McCormick, Davide Ruggero, and Nahum Sonenberg*
October 27-30, 2016 • San Francisco, CA

DNA Repair: Tumor Development and Therapeutic Response

Conference Co-Chairpersons: *Robert G. Bristow, Maria Jasin, and Theodore S. Lawrence*
November 2-5, 2016 • Montreal, Quebec, Canada

New Horizons in Cancer Research: Delivering Cures Through Cancer Science

Conference Co-Chairpersons: *José Baselga and Scott A. Armstrong*
November 2-5, 2016 • Shanghai, P.R. China

Improving Cancer Risk Prediction for Prevention and Early Detection

Conference Co-Chairpersons: *Graham A. Coditz, Susan M. Gapstur, Kenneth R. Muir, and Mark E. Sherman*
November 16-19, 2016 • Orlando, FL

EORTC-NCI-AACR Molecular Targets and Cancer Therapeutics Symposium

Conference Co-Chairpersons: *Jean-Charles Soria, Lee J. Helman, and Levi A. Garraway*
November 29-December 2, 2016 • Munich, Germany

San Antonio Breast Cancer Symposium

Symposium Co-Directors: *Carlos L. Arteaga, Virginia G. Kaklamani, and C. Kent Osborne*
December 6-10, 2016 • San Antonio, TX

Precision Medicine Series: Opportunities and Challenges of Exploiting Synthetic Lethality in Cancer

Conference Co-Chairpersons: *René Bernards, William C. Hahn, and Louis M. Staudt*
January 4-7, 2017 • San Diego, CA

AACR International Conference on New Frontiers in Cancer Research

Conference Co-Chairpersons: *Peter A. Jones and Frank McCormick*
January 18-22, 2017 • Cape Town, South Africa

AACR Annual Meeting 2017

Program Committee Chairperson: *Kornelia Polyak*
April 1-5, 2017 • Washington, DC

Learn more and register at
AACR.org/Calendar

AACR^American Association
for Cancer Research

FINDING CURES TOGETHER™

Announcing Keystone Symposia's 2016–2017 RNA Conferences

Precision Genome Engineering

Scientific Organizers: J. Keith Joung, Emmanuelle Charpentier and Olivier Danos

January 8–12, 2017 | Breckenridge, Colorado | USA

www.keystonesymposia.org/17A2

Deadlines: Scholarship/Discounted Abstract – Sep 13, 2016; Abstract – Oct 6, 2016; Discounted Registration – Nov 9, 2016

Omics Strategies to Study the Proteome

Scientific Organizers: Alan Saghatelian, Chuan He and Ileana M. Cristea

January 29–February 2, 2017 | Breckenridge, Colorado | USA

www.keystonesymposia.org/17A8

Deadlines: Scholarship/Discounted Abstract – Sep 28, 2016; Abstract – Oct 27, 2016; Discounted Registration – Nov 29, 2016

Epigenetics and Human Disease: Progress from Mechanisms to Therapeutics

Scientific Organizers: Johnathan R. Whetstone, Jessica K. Tyler and Rabinder K. Prinjha

January 29–February 2, 2017 | Seattle, Washington | USA

www.keystonesymposia.org/17A9

Deadlines: Scholarship/Discounted Abstract – Sep 29, 2016; Abstract – Oct 27, 2016; Discounted Registration – Nov 30, 2016

Noncoding RNAs from Disease to Targeted Therapeutics

Scientific Organizers: Kevin V. Morris, Archa Fox and Paloma Hoban Giangrande

joint with **Protein-RNA Interactions: Scale, Mechanisms, Structure and Function of Coding and Noncoding RNPs**

Scientific Organizers: Gene W. Yeo, Jernej Ule, Karla Neugebauer and Melissa J. Moore

February 5–9, 2017 | Banff, Alberta | Canada

www.keystonesymposia.org/17J5 | www.keystonesymposia.org/17J6

Deadlines: Scholarship/Discounted Abstract – Oct 5, 2016; Abstract – Nov 2, 2016; Discounted Registration – Dec 6, 2016

mRNA Processing and Human Disease

Scientific Organizers: James L. Manley, Siddhartha Mukherjee and Gideon Dreyfuss

March 5–8, 2017 | Taos, New Mexico | USA

www.keystonesymposia.org/17C3

Deadlines: Scholarship/Discounted Abstract – Nov 2, 2016; Abstract – Dec 6, 2016; Discounted Registration – Jan 10, 2017

RNA-Based Approaches in Cardiovascular Disease

Scientific Organizers: Thomas Thum and Roger J. Hajjar


joint with **Molecular Mechanisms of Heart Development**

Scientific Organizers: Benoit G. Bruneau, Brian L. Black and Margaret E. Buckingham

March 26–30, 2017 | Keystone, Colorado | USA

www.keystonesymposia.org/17X8


Deadlines: Scholarship/Discounted Abstract – Nov 30, 2016; Abstract – Jan 11, 2017; Discounted Registration – Jan 26, 2017



Submit an abstract to participate fully in the conference via a poster presentation and possible selection for a short talk. Scholarships are available for graduate students and postdoctoral fellows. For full program, speaker, abstract and scholarship details, visit www.keystonesymposia.org/genetics

KEYSTONE SYMPOSIA™
on Molecular and Cellular Biology
Accelerating Life Science Discovery

PO Box 1630 | Silverthorne, CO 80498 | 800-253-0685 | 970-262-1230



CRISPR-Cas

A Laboratory Manual



The essential guide to CRISPR-Cas

Edited by Jennifer Doudna, *University of California, Berkeley*;
Prashant Mali, *University of California, San Diego*

The development of CRISPR-Cas technology is revolutionizing biology. Based on machinery bacteria use to target foreign nucleic acids, these powerful techniques allow investigators to edit nucleic acids and modulate gene expression more rapidly and accurately than ever before.

Featuring contributions from leading figures in the CRISPR-Cas field, this laboratory manual presents a state-of-the-art guide to the technology. It includes step-by-step protocols for applying CRISPR-Cas-based techniques in various systems, including yeast, zebrafish, *Drosophila*, mice, and cultured cells (e.g., human pluripotent stem cells). The contributors cover web-based tools and approaches for designing guide RNAs that precisely target genes of interest, methods for preparing and delivering CRISPR-Cas reagents into cells, and ways to screen for cells that harbor the desired genetic changes. Strategies for optimizing CRISPR-Cas in each system—especially for minimizing off-target effects—are also provided.

Authors also describe other applications of the CRISPR-Cas system, including its use for regulating genome activation and repression, and discuss the development of next-generation CRISPR-Cas tools. The book is thus an essential laboratory resource for all cell, molecular, and developmental biologists, as well as biochemists, geneticists, and all who seek to expand their biotechnology toolkits.

2016, 192 pages, illustrated (20 color, 4 B&W), index

Paperback: Print Book + eBook \$210

Print Book \$110

eBook \$100

ISBN 978-1-621821-31-1

Hardcover: Print Book + eBook \$250

Print Book \$150

eBook \$100

ISBN 978-1-621821-30-4

Visit our website for special sale pricing!

eBook available exclusively at www.cshlpress.org

Contents

Preface

CHAPTER 1 Overview of CRISPR-Cas9 Biology

INTRODUCTION

Overview of CRISPR-Cas9 Biology
Hannah K. Ratner, Timothy R. Sampson, and David S. Weiss

CHAPTER 2 Guide RNAs: A Glimpse at the Sequences that Drive CRISPR-Cas Systems

INTRODUCTION

Guide RNAs: A Glimpse at the Sequences that Drive CRISPR-Cas Systems
Alexandra E. Briner and Rudolphe Barrangou

PROTOCOL

Prediction and Validation of Native and Engineered Cas9 Guide Sequences
Alexandra E. Briner, Emily D. Henriksen, and Rudolphe Barrangou

CHAPTER 3 Characterization of Cas9-Guide RNA Orthologs

INTRODUCTION

Characterization of Cas9-Guide RNA Orthologs
Jonathan L. Braff, Stephanie J. Yaung, Kevin M. Esvelt, and George M. Church

PROTOCOL

Characterizing Cas9 Protospacer-Adjacent Motifs with High-Throughput Sequencing of Library Depletion Experiments
Jonathan L. Braff, Stephanie J. Yaung, Kevin M. Esvelt, and George M. Church

CHAPTER 4 Large-Scale Single Guide RNA Library Construction and Use for CRISPR-Cas9-Based Genetic Screens

INTRODUCTION

Large-Scale Single Guide RNA Library Construction and Use for CRISPR-Cas9-Based Genetic Screens
Tim Wang, Eric S. Lander, and David M. Sabatini

PROTOCOLS

Single Guide RNA Library Design and Construction

Tim Wang, Eric S. Lander, and David M. Sabatini

Viral Packaging and Cell Culture for CRISPR-Based Screens

Tim Wang, Eric S. Lander, and David M. Sabatini

CHAPTER 5 Adeno-Associated Virus-Mediated Delivery of CRISPR-Cas Systems for Genome Engineering in Mammalian Cells

PROTOCOL

Adeno-Associated Virus-Mediated Delivery of CRISPR-Cas Systems for Genome Engineering in Mammalian Cells
Thomas Gaj and David V. Schaffer

CHAPTER 6 Detecting Single-Nucleotide Substitutions Induced by Genome Editing

INTRODUCTION

Detecting Single-Nucleotide Substitutions Induced by Genome Editing
Yuichiro Miyaoka, Amanda H. Chan, and Bruce R. Conklin

PROTOCOL

Using Digital Polymerase Chain Reaction to Detect Single-Nucleotide Substitutions Induced by Genome Editing
Yuichiro Miyaoka, Amanda H. Chan, and Bruce R. Conklin

CHAPTER 7 CRISPR-Cas9 Genome Engineering in *Saccharomyces cerevisiae* Cells

PROTOCOL

CRISPR-Cas9 Genome Engineering in *Saccharomyces cerevisiae* Cells
Owen W. Ryan, Snigdha Poddar, and Jamie H.D. Cate

CHAPTER 8 Cas9-Mediated Genome Engineering in *Drosophila melanogaster*

INTRODUCTION

Cas9-Mediated Genome Engineering in *Drosophila melanogaster*

Benjamin E. Housden and Norbert Perrimon

PROTOCOLS

Design and Generation of Donor Constructs for Genome Engineering in *Drosophila*

Benjamin E. Housden and Norbert Perrimon

Detection of Indel Mutations in *Drosophila* by High-Resolution Melt Analysis (HRMA)

Benjamin E. Housden and Norbert Perrimon

Design and Generation of *Drosophila* Single Guide RNA Expression Constructs

Benjamin E. Housden, Yanhui Hu, and Norbert Perrimon

CHAPTER 9 Optimization Strategies for the CRISPR-Cas9 Genome-Editing System

INTRODUCTION

Optimization Strategies for the CRISPR-Cas9 Genome-Editing System

Charles E. Vejnar, Miguel Moreno-Mateos, Daniel Cifuentes, Ariel A. Bazzini, and Antonio J. Giraldez

PROTOCOL

Optimized CRISPR-Cas9 System for Genome Editing in Zebrafish

Charles E. Vejnar, Miguel Moreno-Mateos, Daniel Cifuentes, Ariel A. Bazzini, and Antonio J. Giraldez

CHAPTER 10 Editing the Mouse Genome Using the CRISPR-Cas9 System

INTRODUCTION

Editing the Mouse Genome Using the CRISPR-Cas9 System

Adam Williams, Jorge Henao-Mejia, and Richard A. Flavell

PROTOCOL

Generation of Genetically Modified Mice Using the CRISPR-Cas9 Genome-Editing System

Jorge Henao-Mejia, Adam Williams, Anthony Rongvaux, Judith Stein, Cynthia Hughes, and Richard A. Flavell

CHAPTER 11 Genome Editing in Human Pluripotent Stem Cells

INTRODUCTION

Genome Editing in Human Pluripotent Stem Cells

Cory Smith, Zhaohui Ye, and Linzhao Cheng

PROTOCOL

A Method for Genome Editing in Human Pluripotent Stem Cells

Cory Smith, Zhaohui Ye, and Linzhao Cheng

CHAPTER 12 An Introduction to CRISPR Technology for Genome Activation and Repression in Mammalian Cells

INTRODUCTION

An Introduction to CRISPR Technology for Genome Activation and Repression in Mammalian Cells

Dan Du and Lei S. Qi

PROTOCOL

CRISPR Technology for Genome Activation and Repression in Mammalian Cells

Dan Du and Lei S. Qi

APPENDIX General Safety and Hazardous Material Information

INDEX




www.cshlpress.org



Decoding the Language of Genetics

By David Botstein, Lewis-Sigler Institute for Integrative Genomics

In this book, the distinguished geneticist David Botstein offers help and advice to scientists and physicians daunted by the arcane technical terms that flourish in his discipline. The science of gene function has a vocabulary of specialized, sometimes confusing terms to explain how traits and diseases are inherited, how genes are organized and regulated in the genome, and how the genetic code is read and translated by cells. These terms are often a barrier to full understanding of the underlying concepts. Yet, as more and more individuals learn about their genomes, the information these sequences contain cannot be understood or explained without reference to the basic ideas of genetics. Botstein draws on his long experience as a teacher and pioneering scientist to explain and illuminate what many genetic terms mean and how they entered common usage.



2015, 240 pages, illustrated (30 4C, 10 B&W), index
Hardcover \$79


ISBN 978-1-621820-92-5

Contents

Preface	9 Functional Suppression
Acknowledgments	10 The Genetics of Complex Phenotypes
1 The Basics	11 Transcriptional Regulation of Gene Expression
2 Implicit Experiments and the Functional Gene	12 The Modular Architecture of Genes and Genomes
3 Recombination and Linkage Mapping	13 Evolution Conserves Functional Sequences
4 Pathway Analysis	14 Human Population Genetics
5 Regulation of Metabolic Pathways	15 Inferring Human Gene Function from Disease Alleles
6 Phage and the Beginning of Molecular Genetics	16 What Is Next in Genetics and Genomics?
7 Transcription, Translation, and the Genetic Code	Index
8 Suppression Genetics	



www.cshlpress.org



LEAD NATIONAL-SCALE GENOMICS PROJECTS

We are scientists, engineers, communicators, and data lovers.
We are musicians, homebrewers, yogis, cyclists, and chess players.

Our team is passionate about how genomics and technology can impact and improve human health. And, as part of an innovative startup building the most advanced cloud computing platform for biomedical data analysis - that impact is huge.


**At Seven Bridges, turn your ideas, research,
and software into tools used by researchers
around the world.**

SevenBridges

jobs.sevenbridges.com

Immunoassays for immuno-oncology research

Move confidently across platforms



Immune checkpoints have come to the forefront of cancer therapies as a powerful and promising strategy to stimulate anti-tumor T cell activity. Affymetrix provides a wide range of products for immune-oncology research for the quantification of these important checkpoint proteins. ProcartaPlex® Human Immuno-oncology Checkpoint Panel enables simultaneous quantification of 14 checkpoint modulators. Speed, sensitivity, and specificity make Platinum ELISA the method of choice for biomarker evaluation.

Affymetrix offers multiple immunoassay platforms, providing the freedom to move from a multiplex assay to an ELISA while producing comparable data.

- 98% correlation of antibody pairs across platforms
- Lot-to-lot consistency for reliable reproducibility
- Exceptional scalability between ProcartaPlex® assays and traditional Platinum ELISA

Visit our website to learn more about the ProcartaPlex® Human Immuno-oncology Checkpoint Marker Panel and Platinum ELISA

www.ebioscience.com/immuno-oncology

© 2016 Affymetrix, Inc. All rights reserved. For Research Use Only. Not for use in diagnostic procedures.