



Quantitation of specific transcripts by RT-PCR SNUPE assay.

J Singer-Sam

Genome Res. 1994 3: S48-S50

References This article cites 4 articles, 1 of which can be accessed free at:
<http://genome.cshlp.org/content/3/4/S48.full.html#ref-list-1>

License

Email Alerting Service Receive free email alerts when new articles cite this article - sign up in the box at the top right corner of the article or [click here](#).

A horizontal banner advertisement with a teal background. On the left, the text "CRISPR and RNAi Genetic Screening. Your new superpower." is written in white. In the center, there is a white-bordered box containing the words "LEARN MORE" in blue. On the right, there is a photograph of a woman wearing a red and white superhero costume with a red mask. To her right is the Cellecta logo, which consists of a cluster of green dots of varying sizes, with the word "CELLECTA" in white capital letters below it.

To subscribe to *Genome Research* go to:
<https://genome.cshlp.org/subscriptions>

Copyright © Cold Spring Harbor Laboratory Press

Quantitation of Specific Transcripts by RT-PCR SNuPE Assay

Judith Singer-Sam

Department of Biology, Beckman Research Institute of the City of Hope, Duarte, California 91010

The single nucleotide primer extension assay (SNUPE) provides a means of identifying and quantifying allelic variants in DNA and RNA.^(1,2) It has many uses, including the measurement of specific RNA levels, as well as the study of allelic variation and allele-specific expression. The SNUPE assay is done after amplification by PCR of a segment containing the allelic difference, which may be even a single nucleotide. The assay consists of the enzymatic extension by 1 base of a primer just 5' to the position of mismatch in the presence of a [³²P]dNTP specific for either the wild-type or variant sequence. The ratio of radioactivity in the primer products of these two reactions is then determined after denaturing gel electrophoresis.

A major advantage of the method is its usefulness for quantitative measurement over a wide range. The method allows detection of a given transcript in up to 1000-fold excess of RNA from the other allele, depending on which nucleotides differ. The method can be readily adapted for quantitation of absolute amounts of a specific transcript by addition of an internal standard.⁽³⁾ The protocol described here is the one we use for quantitation of the ratio of *Pgk-1a/Pgk-1b* transcripts in total RNA.

REAGENTS

Reagents are as follows:

1. 10× RT-PCR buffer (100 mM Tris-HCl at pH 8.3, 500 mM KCl). 25 mM MgCl₂.
2. 10× SNUPE buffer (100 mM Tris-HCl at pH 8.3, 500 mM KCl, 20 mM MgCl₂, 0.01% gelatin).⁽²⁾
3. TE (10 mM Tris-HCl, 1 mM EDTA at pH 8.0).
4. RNasin, 40 U/μl (Promega, Madison, WI).
5. Moloney murine leukemia virus reverse transcriptase, 200 U/μl (BRL, Gaithersburg, MD).
6. AmpliTaq polymerase 5 U/μl (Cetus, Norwalk, CT).
7. Gelase, 1 U/μl (Epicentre Technologies, Madison WI).
8. Upstream and downstream primers yielding an amplified RT-PCR product, which includes the region of mismatch (20 μM).
9. SNUPE primer, 18-mer or longer, with the 3' end just 5' to mismatched base (20 μM).⁽²⁾
10. 2.5 mM dNTPs diluted in H₂O from 100 mM stock solutions (Boehringer Mannheim, Indianapolis, IN).
11. [³²P]dNTPs 3000 Ci/mmmole, corresponding to the mismatched bases, diluted to 2 μCi/μl in H₂O just before incubation.
12. 20 mg/ml of mussel glycogen (Sigma, St. Louis, MO).
13. Phenol.
14. Chloroform-isoamyl alcohol (24:1).
15. 7.5 M ammonium acetate.
16. Low melting agarose, ultra pure (BRL).
17. 15% polyacrylamide-7 M urea gel. Ethanol.

PROTOCOL

RT-PCR⁽⁴⁾

1. Mix together at 4°C in a final volume of 20 μl: 1–100 ng of RNA, 1× RT-PCR buffer, 5 mM MgCl₂, 20 units of RNasin, 1 μM downstream primer, 1 mM each dNTP, and 50 units of reverse transcriptase. (For more than one reaction, master mix containing all reagents may be used).
2. Incubate for 15–30 min at 42°C.

3. Heat to 99°C for 5 min.
4. Cool to 50°C for hot start.⁽⁵⁾
5. To each tube add 80 μ l of a prewarmed master mix containing RT-PCR buffer, MgCl₂, upstream primer, and AmpliTaq polymerase, so that the final reaction contains, in a volume of 100 μ l, 10 mM Tris-HCl at pH 8.3, 50 mM KCl, 2 mM MgCl₂, 200 μ M dNTPs, 0.2 μ M upstream and downstream primer and AmpliTaq polymerase, 2.5 units per reaction.
6. Add 100 μ l of mineral oil to each tube.
7. Carry out PCR in a thermal cycler for a sufficient number of cycles to give a strong band on an ethidium bromide-stained gel after amplification.

Purification of AmpliTaq Product

1. Extract RT-PCR products with phenol and chloroform/isoamyl alcohol, and precipitate with ethanol after addition of ½ volume 7.5 M ammonium acetate.⁽⁶⁾ Add glycogen, 20–40 μ g, if needed to obtain a visible pellet.
2. After resuspension of the pellet in TE, run samples on low melting agarose gel, and extract with Gelase, following instructions of the manufacturer.
3. Precipitate with ½ volume of 7.5 M ammonium acetate and 2 volumes of ethanol. Add glycogen, 20–40 μ g, as a carrier if needed.
4. Resuspend the DNA in 10–20 μ l of TE, and estimate the concentration of an aliquot on an ethidium bromide-stained gel.

Quantitative SNUPE Assay^(1,2)

1. To each of two tubes, add 10 ng of template DNA, 1 \times SNUPE buffer, 1 μ M SNUPE primer, 0.75 units of *Taq* polymerase, and 2 μ Ci of the [³²P]dNTP specific for one allele or the other (final volume, 10 μ l).
2. Incubate the samples for one round of denaturation, annealing and synthesis (example: 95°C for 1 min, 42°C for 2 min, 72°C for 1 min).
3. As controls for background and maximum incorporation, incubate 10 ng of amplified product from each allele with the “incorrect” and appropriate [³²P]dNTP, respectively.
4. Electrophorese the samples on a 15% polyacrylamide–urea gel.⁽⁶⁾
5. Determine the amount of radioactivity in each ($n + 1$) band (n = length of SNUPE primer). Occasionally, the base following the incorporated [³²P]dNTP is identical to it, and one sees ($n + 2$) products. In that case, the contribution of both ($n + 1$) and ($n + 2$) products should be included.
6. After subtraction of background, the signal obtained with the [³²P]dNTP is then proportional to the relative amount of each allelic transcript. The background varies with different mismatches. We have found measurement of C in an excess of A (C/A) to be sensitive down to the 0.001% level; other mismatches giving a low background level include C/T, T/C, and A/G. G/C, C/G, and G/A mismatches give a background of ~1%.

ACKNOWLEDGMENTS

I thank Drs. Jeanne M. LeBon and Arthur D. Riggs for their help, and Dr. Pirooska Szabó for kindly sharing unpublished data.

REFERENCES

1. Singer-Sam, J., J.M. LeBon, A. Dai, and A.D. Riggs. 1992. A sensitive, quantitative assay for measurement of allele-specific transcripts differing by a single nucleotide. *PCR Methods Applic.* **1**: 160–163.
2. Kuppaswamy, M.N., J.W. Hoffman, C.K. Kasper, S.G. Spitzer, S.L. Groce, and S.P. Bajaj. 1991.

Manual Supplement

- Single nucleotide primer extension to detect genetic diseases: Experimental application to hemophilia B (factor IX) and cystic fibrosis genes. *Proc. Natl. Acad. Sci.* **88**: 1143–1147.
3. Singer-Sam, J. and A.D. Riggs. 1993. Quantitative analysis of messenger RNA levels: Reverse transcription-polymerase chain reaction single nucleotide primer extension assay. *Methods Enzymol.* **225**: 344–351.
 4. Kawasaki, E.S. 1990. Amplification of RNA. In *PCR protocols*. (ed. M.A. Innis, D.H. Gelfand, J.J. Sninsky, and T.J. White), pp.21–27. Academic Press, San Diego, CA.
 5. D'Aquila, R.T., L.J. Bechtel, J.A. Videler, J.J. Eron, P. Gorczyca, and J.C. Kaplan. 1991. Maximizing sensitivity and specificity of PCR by pre-amplification heating. *Nucleic Acids Res.* **19**: 3749.
 6. Maniatis, T., E.F. Fritsch, and J. Sambrook. 1982. *Molecular cloning: A laboratory manual*. Cold Spring Harbor Laboratory, Cold Spring Harbor, New York.