



**Supplemental Figure S14. Calibration experiments to determine mCherry gating for single integration sorting.** (A) pZDonor plasmid containing a *ACTB* promoter upstream of the eGFP was integrated into K562 using ZFN-mediated site-specific integration. The populations of single integrations (P2) and double integrations (P3) were evaluated according to the GFP and mCherry density plot. (B) Cells were sorted by their mCherry levels, according to differential gating on mCherry (P1(7% of the mCherry positive cells), P5(12.1%), P6(19.3%), P7(26%), P8(32.7%), P9 (40.5%), P10(59.3%), P11 (78.2%). As a control, the range containing both single integrations and double integrations was sorted as well (P4(7.4%)). (C) Cells were reanalyzed 7 days after the sorting according to mCherry levels. While all the gates up to P9 (P1,P5,P9 are shown) exhibited only 1% or less of the cells in the P3 gate, i.e. cells with double integration, in P10 1.8% of the cells were present in this gate, and 9.5% in P11. P9 (highlighted in a red frame) was thus chosen as the gate allowing collection of most cells, with minimal rate of cells that may possess double integrations. In the positive control P4, 36% of the cells were in P3 a week after the sorting.